InVivoMAb Armenian hamster IgG isotype control, anti-glutathione S-tr

Lot Specific Information

Lot Number:	Lot Specific*
Volume:	Lot Specific*
Concentration:	Lot Specific* (generally 4 to 11 mg/ml) *
Total Protein:	Lot Specific*
*This information will be noted on the certificate of analysis that ships with this product.	

Product Information

Catalog Number:	BE0260
Clone:	PIP
Isotype:	Armenian hamster IgG
Recommended Dilution Buffer:	InVivoPure pH 7.0 Dilution Buffer
Immunogen:	
Reported Applications:	
Formulation:	PBS, pH 7.0 Contains no stabilizers or preservatives
Endotoxin:	<2EU/mg (<0.002EU/µg) Determined by LAL gel clotting assay
Purity:	>95% Determined by SDS-PAGE
Sterility:	0.2 µM filtered
Production:	Purified from tissue culture supernatant in an animal free facility
Purification:	Protein A
RRID:	AB_2687739
Molecular Weight:	150 kDa

Description

The PIP monoclonal antibody reacts with glutathione-S-transferase (GST) of *Schistosoma japonicum* origin. Because *S. japonicum* GST is not expressed by mammals this antibody is ideal for use as an isotype control for Armenian hamster IgG antibodies in most *in vivo* and *in vitro* applications.

Shelf-life and Storage

Store at the stock concentration at 4°C. Do not freeze.

All Bio X Cell antibodies have a guaranteed shelf-life of one year from the date of customer receipt when stored as recommended. It is not uncommon for a floccule or precipitate to appear during storage. The floccule is typically buffer salts precipitating out of solution or a small bit of protein aggregation. For information on how to remove floccules or precipitates see our FAQ's at <u>bxcell.com/faqs</u>.

Protocol Information

Since applications vary, each investigator should use the application references as a guide to help estimate the appropriate dose or concentration. The dose or concentration can be further optimized experimentally in a dose response or titration experiment.

Application References

For a complete list of references, visit https://bxcell.com/product/invivomab-anti-qst/#references_or scan the QR code below.



Bio X Cell, Inc.

Bio X Cell, Inc.