



**InVivoMAb human IgG2 isotype control**

**Lot Specific Information**

<b>Lot Number:</b>	Lot Specific*
<b>Volume:</b>	Lot Specific*
<b>Concentration:</b>	Lot Specific* (generally 4 to 11 mg/ml) *
<b>Total Protein:</b>	Lot Specific*

\*This information will be noted on the certificate of analysis that ships with this product.

**Product Information**

<b>Catalog Number:</b>	<b>BE0301</b>
<b>Clone:</b>	<b>N/A</b>
<b>Isotype:</b>	Human IgG2, λ
<b>Recommended Dilution Buffer:</b>	InVivoPure pH 7.0 Dilution Buffer
<b>Immunogen:</b>	
<b>Reported Applications:</b>	
<b>Formulation:</b>	PBS, pH 7.0 Contains no stabilizers or preservatives
<b>Endotoxin:</b>	<2EU/mg (<0.002EU/μg) Determined by LAL gel clotting assay
<b>Purity:</b>	>95% Determined by SDS-PAGE
<b>Sterility:</b>	0.2 μm filtration
<b>Production:</b>	Purified from human myeloma serum
<b>Purification:</b>	Protein A
<b>RRID:</b>	AB_2715459
<b>Molecular Weight:</b>	150 kDa

**Human Pathogen Test Results**

Human immunodeficiency virus: Negative, Hepatitis B virus: Negative, Hepatitis C virus: Negative,

**Description**

The human IgG2 isotype control antibody is purified from human myeloma serum and is of unknown specificity. This antibody is suitable for use as a non-targeting isotype control in various *in vitro* and *in vivo* studies. It can also be used as a negative control in various diagnostic applications such as ELISA, Western blot, immunofluorescence, immunohistochemistry, immunoprecipitation, and flow cytometry. For research use only.

**Shelf-life and Storage**

Store at the stock concentration at 4°C. **Do not freeze.**  
All Bio X Cell antibodies have a guaranteed shelf-life of one year from the date of customer receipt when stored as recommended. It is not uncommon for a floccule or precipitate to appear during storage. The floccule is typically buffer salts precipitating out of solution or a small bit of protein aggregation. For information on how to remove floccules or precipitates see our FAQ's at [bxcell.com/faqs](https://bxcell.com/faqs).

**Protocol Information**

Since applications vary, each investigator should use the application references as a guide to help estimate the appropriate dose or concentration. The dose or concentration can be further optimized experimentally in a dose response or titration experiment.

**Application References**

For a complete list of references, visit <https://bxcell.com/product/invivomab-human-igg2-isotype-control/#references> or scan the QR code below.

**Bio X Cell, Inc.**

bxcell.com  
1.866.787.3444  
[customerservice@bxcell.com](mailto:customerservice@bxcell.com)

Conditions: For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.  
Bio X Cell, Bio X Cell Logo and all other trademarks are the property of Bio X Cell, Inc. © 2020 Bio X Cell



**Bio X Cell, Inc.**

---

[bxcell.com](http://bxcell.com)

1.866.787.3444

[customerservice@bxcell.com](mailto:customerservice@bxcell.com)

*Conditions: For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.*

**Bio X Cell, Bio X Cell Logo and all other trademarks are the property of Bio X Cell, Inc. © 2020 Bio X Cell**